

Tracing of visual projections

RV Robin J Vigouroux AC Alain Chédotal KN Kim Tuyen Nguyen-Ba-Charvet

Updated date: Jan 16, 2023



An abbreviated version of this protocol was published in eLIFE in Feb 2020

Revisiting the role of Dcc in visual system development with a novel eye clearing method

DOI: 10.7554/eLife.51275

Detailed protocol

Dear Carlos,

For any species (mice, fish), we used AlexaFluor555-CTB C22843 or AlexaFluor647-CTB C34778 at 0.2 % diluted in water. We then aliquot in 0.5 ml microtubes and store them at -20°C. As we're using a lot of CTB, we never keep them at -20°C more than a few months.

Kim

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Vigouroux, R. , Chédotal, A. and Nguyen-Ba-Charvet, K. (2023). Tracing of visual projections. Bio-protocol Preprint. bio-protocol.org/prep2119.
2. Vigouroux, R. J., Cesar, Q., Chédotal, A. and Nguyen-Ba-Charvet, K. T.(2020). Revisiting the role of Dcc in visual system development with a novel eye clearing method. eLIFE. DOI: [10.7554/eLife.51275](https://doi.org/10.7554/eLife.51275)

Copyright: Content may be subjected to copyright.